NEW ENDOPREDICT RESULT REPORT QUICK GUIDE



One Test - Three Clinical Answers

EndoPredict is the only test that answers three critical clinical questions:


- 1. Risk at 10 years Can chemotherapy be avoided?
- 2. Chemotherapy benefit What is the absolute benefit from chemotherapy?
- 3. Risk between 5 and 15 years Can the extension of endocrine therapy be avoided?

EndoPredict is a second-generation gene expression assay that combines a molecular score with clinical features (tumor size and nodal status). The five-page report provides the risk of developing distant recurrence (0-10 and 5-15 years) and the absolute chemotherapy benefit at 10 years for a woman with ER+, HER2- breast cancer. Results are summarized by treatment planning stage on the front page. Detailed information is available on the pages that follow.

The EndoPredict Result Report uses the colors of the arrows of the product logo to provide a clear and immediate indication of the patient's risk: blue EPclin Risk Score and EPclin Risk Class indicate high risk, while green indicates low risk.

Elements of the EndoPredict Result Report

Front Page



NEW ENDOPREDICT RESULT REPORT QUICK GUIDE




Page 2

LIKELIHOOD OF DISTANT RECURRENCE WITHIN YEARS 0-10

EPclin Risk Score's continuous risk curve

 Graphically illustrates the patient's individualized risk of distant recurrence within years 0-10 based on the ABCSG 6&8 validation studies.



10-year likelihood of distant recurrence for the individual patient

 Patient's individual risk of distant recurrence when treated with 5 years of endocrine therapy alone.

Since launch, EndoPredict has provided an accurate answer to an important clinical question: What is the risk of distant recurrence in the first ten years after diagnosis with endocrine therapy only? With this information, your patients with a low risk of recurrence could be identified and could safely avoid chemotherapy.

Page 4

LIKELIHOOD OF LATE DISTANT RECURRENCE WITHIN YEARS 5-15

EPclin Risk Score's continuous risk curve

 Graphically illustrates the patient's individualized risk of distant recurrence within years 5-15 based on the ABCSG 6&8 validation studies.



10-year likelihood of late distant recurrence for the individual patient

 Patient's individual late distant recurrence risk. Provides estimates of distant recurrence at 15 years for patients who are distant recurrence-free after 5 years of endocrine therapy alone and who do not receive extended endocrine therapy.


The third answered question is: What is your patient's risk of late recurrence (5-15 years)? This information guides treatment decisions to extend endocrine therapy beyond 5 years.

Page 3

ABSOLUTE CHEMOTHERAPY BENEFIT AT 10 YEARS

EPclin Risk Score's continuous risk curves to determine the — chemotherapy benefit

 Graphically illustrates the patient's individualized chemotherapy benefit based on the ABCSG 6&8, TransATAC, GEICAM/2003-02 and GEICAM/9906 studies.



Absolute chemotherapy benefit for the individual patient

· Patient's individual absolute chemotherapy benefit.

With the new report EndoPredict answers a second question: What is the personalized absolute chemotherapy benefit? This leads to a confident chemotherapy treatment decision for your patient.

Page 5

Page 5 of the **Result Report** provides additional details on the EndoPredict assay design, intended use, and its clinical applications.

Important Note

An EndoPredict Result Report is only generated if both patient's tumor size and nodal status are available. If an EndoPredict test is performed on a biopsy specimen, and therefore patient's tumor pathological information is not yet available, the EndoPredict Result Report is not provided. However, the 12-Gene Molecular Score and a table showing the EPclin Risk Score and the EPclin Risk Class for all possible combinations of Tumor Size and Nodal Status will be available in the EndoPredict Table Report. If you would like an updated report with all risk curves you can send us the tumor size and nodal status after resection and request an amended report.

For additional information regarding the EndoPredict Result and Summary Reports, please refer to the EndoPredict Quick Guides and the EndoPredict Instruction Manuals, available from Myriad International. Please contact: epsupport@myriadgenetics.eu.